

尿液培養與敏感性試驗檢查

您的醫務人員要求您提供一個尿液的樣本，以進行尿液培養與敏感性試驗檢查的測試。為了獲得有意義的實驗室測試結果，請嚴格遵循這些說明並請從我們的實驗室獲取無菌樣品容器。

採集尿液

在起床後或任何時間，將尿液收集到所提供的 C&S 容器內。

1. 洗淨並擦乾雙手。
2. 從包裝中取出濕巾，用濕巾一次性從前到後清潔陰莖或陰唇。
3. 開始向馬桶小便少量。把容器放在尿流處收集尿液時，不要停止尿流，在尿流結束前移開容器。

第 3 步 中段尿圖解：



4. 緊緊關閉容器蓋，注意您的手指不要觸摸到容器的邊緣和內側。
5. 用肥皂和水洗手。

貼標籤

1. 在容器上所貼的標籤要有您的姓名、出生日期以及採樣日期。

包裝

1. 將盛尿的容器放入提供的袋子，並牢固地封閉袋口。
2. 把所有書面文件疊好，放入樣本運輸袋的外置兜內。這就確保樣本不會泄漏到文件上。

儲存與運輸

1. 把樣本冷藏，並在採樣的當日把樣本送到實驗室。

Revision Number:	1.0	Effective Date:	Nov 16, 2022
------------------	-----	-----------------	--------------

NOTE: This is a **CONTROLLED** document. Any documents appearing in paper form that are not stamped in red "MASTER COPY" are not controlled.

PRE-PR-0024

Prepared by: Bonnie Tam

Approved By:

Medical Approval		
Dr. Renison Chongkit		Nov 16, 2022
Name	Signature	Date

Management Approval		
Bonnie Tam		Nov 16, 2022
Name	Signature	Date

Record of Revisions

Version	Date of Revision/Review	Responsible Person	Revision Description	Signature of Approval

Revision Number:	1.0	Effective Date:	Nov 16, 2022
-------------------------	-----	------------------------	--------------

*NOTE: This is a **CONTROLLED** document. Any documents appearing in paper form that are not stamped in red "MASTER COPY" are not controlled.*